





A new generation light weight aluminium lifter evolved from the needs of the carer.



MICROS | 30







Quick release mast





Simple adjustment of base width with foot pedal





Extra low 34cm lifting height to allow floor lifts





B

PRODUCT DESCRIPTION

Compact Aluminium lift Micros 7108.

- Aluminium lift with a Safe Working Load (SWL) of 130kg
- Linak powered lifting function •
- Foot pedal control for base width adjustment
- Strong and comfortable handles for easy manoeuvring
- Equipped with 75mm front castors and 75mm lockable rear castors •
- Jumbo Care diagnostic control and service monitor system
- Easy access emergency stop button •
- Electric and manual emergency lowering devices •
- Anti-entrapment feature lowering of boom stops on resistance





MAX MIN MAX MIN						-	U	L	B1	B	4		1		MAX LOAD	ITEM NO
											MIN	MAX	MIN	MAX		
7108 130 151 30 173 103 104 82 57.6~94.2 47.4~84.2 94.5 11.2 5.4 89.5 31kg	117	31kg	89.5	5.4	11.2	94.5	47.4~84.2	57.6~94.2	82	104	103	173	30	151	130	7108

Il measurements are shown in centimeth

Allegro Concepts Pty Ltd

30-32 James Street Lidcombe NSW 2141 Phone: (02) 9749 7812 Fax (02) 9749 2144 Email: sales@allegroconcepts.com.au

