

Hoist Attachments - General Purpose



Standard Yoke

- 2 sling attachment points at extremity of curved yoke bar
- Use with general purpose (loop attachment) slings
- Non-absorbent polyurethane padding

Available on: Alto, Sonata & Sonata ELS

Indications: Suitable for the majority of patients

Useful for: All day-to-day patient lifting



Standard Yoke - SY02

4 Point Yoke

- 4 sling attachment points spread over a broad area
- Use with general purpose (loop attachment) slings
- Padded vinyl covers with velcro closure

Available on: Alto, Sonata & Sonata ELS

Indications: Bariatric clients, clients requiring a semi-supine lift to maintain upper body control

Useful for: Supine lifting devices eg. Jordan frame



Power Pivot Frame - PPF01

Allegro
CONCEPTS PTY LTD

Innovation in Patient Handling

Hoist Attachments - General Purpose

The Allegro Concepts range of standard patient lifting hoists includes:

- Alto Lift (SWL 200kg)
- Sonata (SWL 150kg)
- Sonata ELS (SWL 150kg)

These three heavy duty hoists are available with a choice of 6 different sling attachment options.

Achieve quality patient transfers by selecting the fitting which best suits your client or healthcare facility.

Factors which may determine choice of sling attachment include:

- Physical attributes of patients eg. weight, disability
- Level of staff training eg. competence with complex lifting procedures
- Necessity to interface with existing equipment eg. beds, wheelchairs

Tri-Hook Yoke

- 6 Sling attachment points at extremity of curved yoke bar
- Use with general purpose (loop attachment) slings
- Non-absorbent polyurethane padding

Available on: Alto, Sonata & Sonata ELS

Indications: Where the sling prescribed has more than the usual 4 attachment loops

Useful for: Facilities which have historically trained with this configuration



Tri Hook Yoke - SY03

Hi-Lift Yoke

- 2 sling attachment points at extremity of slightly curved yoke bar
- Use with general purpose (loop attachment) slings
- Padded vinyl cover with velcro closure

Available on: Alto, Sonata & Sonata ELS

Indications: Where lifting height of specified hoist is insufficient (increases lifting height by 80mm)


Useful for: Reducing the loss of lifting height when installing an Allegro Universal Weighing Device



Hi-Lift Yoke - SY02HL

Bariatric 4 Point Yoke - BY04


- All day-to-day patient lifting
- Easily adaptable to any style of patient hoist
- Optimum patient space for larger patients
- 4 Sling attachment points
- Use with general purpose (loop attachment) slings
- Non-absorbent polyurethane padding
- 320kg SWL



Bariatric 4 Point Yoke - BY04

Bariatric Yoke Spreader Bar with Integrated Weigh Scale

- All day-to-day patient lifting
- Monitor patient wellbeing with no loss of hoist lifting height
- Easily adaptable to any style of patient hoist
- Optimum patient space for larger patients
- 4 Sling attachment points
- Use with general purpose (loop attachment) slings
- Non-absorbent polyurethane padding



Bariatric Yoke Spreader Bar with Integrated Weigh Scale - BY04WD

Allegro Concepts Pty Ltd

Unit 9, 110- 120 Silverwater Road | Silverwater | NSW | 2128

☎ 02 6895 2888 | ✉ sales@allegroconcepts.com.au

🌐 www.allegroconcepts.com.au