

Wireless Cushion Floor Sensor Mat

PREMIUM OUT OF BED CUSHION MAT

THICK FIRM FOAM TO REDUCE, IF NOT ELIMINATE INJURIES

INCLUDE A 2-PIECE SENSOR MAT AND A WIRELESS TRANSMITTER

FOLDABLE FOR EASY STORAGE

EASY TO MOVE FROM ROOM TO ROOM

SPRAY AND WIPE TO CLEAN WITH ANTI-BACTERIA SOLVENT

INTEGRATES TO MOST NURSE CALL SYSTEMS

PLUG & PLAY OPERATION

COST EFFECTIVE

LARGE PRESSURE MAT (1800mm X 750mm X 40mm)

TGA Registered ID 223902
RCM Approved
AS/NZS 61000.6.3.2012



PRODUCT CODE:

Cushion Floor Mat & Transmitter: MED-FM-04X-WT

FULL KIT: Receiver Box, Floor Mat & Transmitter: MED-FM-04K-WT

JB Medical distributes a smart Australian designed Foldable Premium WIRELESS Cushion Floor Mat for residents who are at risk of falling out of their bed. The cordless floor mat purposely designed for easier storage and trip-hazard free. The transmitter is an exclusive ID encoded RF Transmitter which communicates with the Wireless Receiver Box (RRU) or the Universal Call Point (UCPV5.1). It allows the Floor mat to wirelessly transmit an alarm to the call point. Simply connect the Transmitter to the floor sensor mat. Then teach in the Transmitter to RRU or Dual Function Call Point by following the instructions printed. The sensor mat will now trigger the call point when activated within the transmitting area. (Note: Transmitter and Floor Mat can be ordered separately)

KEY SPECIFICATIONS:

- Dimensions: 1800mm X 750mm X 40mm (Foldable)
- Electrical Resistance < 150 Ohm
- Blue PVC Soft Foam Design
- 40 mm of Foam Padding
- Connection: 6.3mm Jack socket
- No External Antenna
- LED low battery indicator
- Transmits up to 30m (line of sight)
- Frequency: 433.92MHz
- Uses a 3V battery-CR2032
- Warranty: 12 months

PRODUCT ACCESSORIES:

- Wireless Receiver (RRU)
- Wireless Bed Mat
- Wireless Under Mattress Mat
- Wireless Chair Mat
- Wireless Waterproof Neck or Wrist Pendant
- Wireless PIR Motion Sensor

JB MEDICAL EQUIPMENT

Unit 1, 4 Progress Circuit, Prestons, NSW, 2170
PH +61 1300-788-182 FAX +61 2-8783-5522
www.jbmedical.com.au
ABN 18114865520



ORDER INFORMATION

Fax or email your purchase orders to
02-8783-5522
sales@jbmedical.com.au

