

Communal Wireless Receiver Unit (RRU)

REDUCE THE RISK OF
FALLS & INJURIES IN
COMMUNAL AREA'S

TURN HARDWIRED
CALL POINTS TO
WIRELESS

UNLIMITED WIRELESS
DEVICES CAN BE USED
IN THAT AREA

NO WIRES AROUND
RESIDENT/PATIENT-NO
TRIP HAZARDS

CHOOSE THE WIRELESS
DEVICES TO SUIT

PORTABLE

EASY TO USE

EASY TO MOVE FROM
ROOM TO ROOM

INTEGRATES TO MOST
NURSE CALL SYSTEMS

PLUG & PLAY
OPERATION

DURABLE & RELIABLE

COST EFFECTIVE



RCM Approved
AS/NZS 61000.6.3.2012

PRODUCT CODE: MED-NC-RRU(COMM)-T(X)-V1.2

X= 1 for Tip (common) and Ring wiring

X= 3 for Tip and Sleeve (common) wiring

Indigo Care's Wireless Communal Receiver Box is a RF Receiver for interfacing to Nurse Call Points in Hospitals and Aged Care Facilities that conform to AS3811.

This device allows unlimited wireless nurse call or fall prevention devices to be used in communal areas.

Wireless devices that are taught into resident's room RRU's, such as wireless pendants on lanyards and chair exit pads on wheel chairs will now also work when the resident moves into the communal area.

Simply plug the Communal RRU into Nurse Call Point socket, there is no need to teach in any devices. Any of our wireless devices will now trigger the call point when activated within the transmitting area.

KEY SPECIFICATIONS:

- Box Dimensions: 67x 67x 22 mm
- Connection: 6.3mm right angle jack plug on 15cm cable
- External Antenna
- Status LED
- Receives up to 25m (line of sight)
- Receives unlimited transmitter ID's
- Operates with a 9VDC Switch mode power pack
- Frequency: 433.92MHz
- Call & Bed light Outputs are Relay switched
- Additional 6.3 Jack and RJ15 outlets for hardwired devices (eg PIR, pendants, floor sensor)
- Warranty: 12 months

PRODUCT ACCESSORIES:

- Wireless Floor Sensor Mat
- Wireless Bed Mat
- Wireless Chair Mat
- Wireless Waterproof pendant
- Wireless PIR Motion Sensor

Horizon Electronic Components

T/as INDIGO CARE

2/71 Rushdale Street, Knoxfield, VIC, 3180

PH +61 3 9763-2988 FAX +61 3 9763-3811

www.indigocare.com.au

ABN 36007302801



ORDER INFORMATION

Fax or email your purchase orders to

03-9763-3811

Or email sales@indigocare.com.au

