

Wireless Bed Sensor Mat

**MORE COMFORT
LESS NOISE**

**INTEGRATES TO MOST
NURSE CALL SYSTEMS
VIA OUR RECEIVERS**

**NO WIRES AROUND
RESIDENT/PATIENT-
NO TRIP HAZARDS**

**REDUCE THE RISK OF
FALLS & INJURIES**

**SOLVENT WIPE-
DISINFECTANT**

**LOW PROFILE, NON-
SLIP MATERIAL**

PORTABLE

**EASY TO MOVE FROM
ROOM TO ROOM**

**PLUG & PLAY
OPERATION**

DURABLE & RELIABLE

COST EFFECTIVE



(TGA Registered ID 223902

RCM Approved

AS/NZS61000.6.3.2012)

PRODUCT CODES: Bed Mat & Transmitter: MED-B365X-WT

Full Kit: Receiver Unit, Bed Mat & Transmitter: MED-B365K-W6T(X)orW3T(X)

JB Medical distributes a New & Improved pressure sensitive Premium Bed Sensor Mat that is more comfortable & less noisy for the user. The mat's thickness has increased to 12mm providing a softer touch and feel for the user at the same time ensuring the noise created due to contact is quieter. This improved Bed Mat will be a huge hit with your residents and staff! MORE COMFORT, LESS NOISE

Our Bed Mats are designed to immediately alert staff should a resident leave their bed which is most useful to residents who are at risk of falling out of their bed, or those who require assistance in arising from the bed. The transmitter is an exclusive ID encoded RF Transmitter which communicates with the Wireless Receiver Box (RRU). It allows the bed sensor mat to wirelessly transmit an alarm to the call point. Simply connect the Transmitter to the bed sensor mat, then teach in the Transmitter to RRU by following the instructions printed. The sensor mat will now trigger the call point when activated within the transmitting area. It can be integrated to most nurse call systems. (Note: Transmitter and Bed mat can be ordered separately)

KEY SPECIFICATIONS:

Bed mat

- **Dimensions: 770mm X 495mm**
- **Connection: RJ11 on 500mm curly cord**
- **Electrical resistance < 500 Ohm**
- **White PVC, Thickness 12mm**
- **Warranty: 12 months**

Transmitter

- **Connection: RJ12 socket**
- **No External Antenna**
- **LED battery indicator**
- **Transmits up to 30m (line of sight)**
- **Frequency: 433.92MHz**
- **Uses a 3V battery-CR2032**
- **Warranty: 12 months**

PRODUCT ACCESSORIES:

- **Wireless Receiver Box**
- **Wireless Floor Sensor Mat**
- **Wireless Chair Mat**
- **Wireless Waterproof pendant**
- **Wireless PIR Motion Sensor**

ORDER INFORMATION

**Fax or email your purchase orders to
02-8783-5522
sales@jbmedical.com.au**

**JB MEDICAL EQUIPMENT
Unit 1, 4 Progress Circuit, Prestons, NSW, 2170
PH +61 1300-788-182 FAX +61 2-8783-5522
www.jbmedical.com.au
ABN 18114865520**

