

Wireless Wrist Pendant

COMBINATION OF A
WRIST & NECK
PENDANT

REDUCE THE RISK OF
FALLS & INJURIES

ALLOWS THE
RESIDENT TO CALL
FOR ASSISTANCE
ANYWHERE WITHIN
THEIR
ROOM/ENSUITE

WATER RESISTANT
IP67

SILICON STRAP BAND

SILICON BUTTON
WITH BRAILLE

RE-ASSURANCE LIGHT

LOW BATTERY
WARNING

INTEGRATES TO MOST
NURSE CALL SYSTEMS
VIA OUR RECEIVERS

DURABLE & RELIABLE

COST EFFECTIVE



PRODUCT CODE: MED-RTU-WT

JB Medical distributes an exclusive Wireless Wrist Pendant which is a combination of a wrist & neck pendant; two assistive technological products in one! Choose whether your resident wears it as a traditional wrist band or as a neck pendant, simply by removing the pendant from the wrist band and attaching the lanyard. This device is an excellent communication tool for times of distress or danger. The right choice of pendant, whether wrist or neck, can give both you and your residents more freedom and security. Easy to use and visually appealing. The pendant is water resistant to IP67 rating. The easy to press CALL button with braille is made of soft silicon. The pendant has a re-assurance LED light for nurse call and low battery warning to the RRU receiver.

The pendant has an exclusive ID encoded RF Transmitter which communicates with the Wireless Receiver Box (RRU) or the Call Point. It enables the resident to call for assistance anywhere in their room or ensuite without the need to move from their resting position. Simply program the Wrist Pendant to the RRU or Call Point by following the printed instructions or view our online tutorials.

The Kit includes: 1 Wristband Pendant, 1 Silicon Neck Band, 2 batteries, 1 snap-joint lanyard

Note: Wrist Pendant is not suitable for use in a spa, sauna, steam or any hydrotherapy area

KEY SPECIFICATIONS:

- **Weight of Neck pendant: 22 grams**
- **Weight of Wrist pendant: 32 grams**
- **Re-assurance Light**
- **Not suitable for use in a spa, sauna, steam or any hydrotherapy area**
- **Water Resistant - IP67**
- **Silicon CALL button with Braille**
- **Low Battery warning to RRU**
- **Transmits up to 30m (line of sight)**
- **Frequency: 433.92MHz**
- **Uses a 3V battery-CR2032**



PRODUCT ACCESSORIES:

- **Wireless Receiver Box - RRU**
- **Wireless Bed Mat**
- **Wireless Chair Mat**
- **Wireless Transmitter-Bed/Chair**
- **Wireless Floor Mat**
- **Wireless PIR Motion Sensor**
- **Wireless Personal Pendant**

JB MEDICAL EQUIPMENT

Unit 1, 4 Progress Circuit, Prestons, NSW, 2170

PH +61 1300-788-182 FAX +61 2-8783-5522

www.jbmedical.com.au

ABN 18114865520

ORDER INFORMATION

Fax or email your purchase orders to

02-8783-5522

Or email sales@jbmedical.com.au

