

# ALTO Lift™ MkII 200

Heavy Duty Patient Lifting Hoist



*Allegro*  
CONCEPTS PTY LTD

Innovation in Patient Handling

# ALTO Lift™ MkII 200

ALTO Lift™ MkII 200: A medium sized heavy duty mobile patient lifting hoist with a safe working load of 200kg and a compact ergonomic design.

ALTO Lift™ MkII 200, whilst being compact, provides superior patient lift space through the use of SPACE LINK™ technology.

ALTO Lift™ MkII 200 is the lifter of choice for institutional use or home care where a high capacity lift is required.

Also available Quick Release Adapter option.



## SPACE LINK™

- breakthrough for optimized patient space
- enhances floor lift capacity

## Alto Lift 200 Specifications

Safe Working Load	200kg
Base Length	1150mm
External Base Width Closed	700mm
External Base Width Open	1130mm
Internal Base Width Closed	600mm
Internal Base Width Open	1030mm
Base Width Adjustment	Electric
Under Bed Clearance	130mm
Minimum Height Yoke	300mm
Maximum Height Yoke	1920mm
Average Lift Range*	1620mm
Lifter Weight	46kg
Minimum Storage Height	1350mm
Under Base Height	55mm
Maximum Hoist Reach	1060mm
Minimum Hoist Reach	610mm
Hoist Reach at Maximum Height	515mm
Mass of Heaviest Part	24kg (Mast assembly)

\* With general purpose large sling

The Allegro Concepts range of patient lifting equipment including lifters, slings and accessories has been designed and tested to comply with: AS ISO 10535-2011 Hoists for the transfer of disabled persons – Requirements and test methods.

## Allegro Concepts Pty Ltd

Unit 9, 110- 120 Silverwater Road | Silverwater | NSW | 2128

☎ 02 6895 2888 | ✉ sales@allegroconcepts.com.au

🌐 www.allegroconcepts.com.au



Standard Yoke - SY02



4 Point Yoke - SY04



Standard Yoke with Integrated Weigh Scale - SY02WD



Pivot Frame - PF01



Powered Pivot Frame - PPF01



Self Power Pivot Frame - SPPF01