



www.jbmedical.com.au

# AUGUST 2022

## BROCHURE

### P5500 BED

**SKU: PHYNIGHTSPLG-5**

**SINGLE: BD0310**

**KING SINGLE: BD0330**

**AVAILABLE IN  
SINGLE & KING  
SINGLE BED**

### SPECIFICATION

Code	Size	Type	Overall Width	Overall Length	Sleeping Surface	Height	Safe Working Load
BD0310	Single	P5500	930mm	2190mm	900 x 2000mm	295-740mm	200kg
BD0330	King Single	P5500W	1070mm	2190mm	1050 x 2000mm	295-740mm	200kg



**PIN SLOT SYSTEM**



**TRANSPORT MODE (FOLDED)**



**HANDSET**

**ON SALE**  
ORDER NOW

### KEY FEATURES

- ✓ ELECTRIC ADJUSTABLE BED IDEAL FOR THE HOMECARE ENVIRONMENT
- ✓ ELECTRIC HIGH/LOW FUNCTIONS, INCLUDING BACKREST AND KNEE-BEND
- ✓ ELECTRIC TRENDELENBURG AND REVERSE TRENDELENBURG POSITIONS
- ✓ PULL-APART FUNCTION FOR EASY TRANSPORT AND INSTALLATION
- ✓ SUPPLIED WITH EUROPEAN BEECH HEAD AND FOOTBOARD PANELS
- ✓ HAND CONTROL IS SIMPLE AND EASY TO USE
- ✓ 4 X 100MM LOCKING CASTORS
- ✓ BEIGE POWDER COATED STEEL FRAME

### ACCESSORIES AVAILABLE

- ✓ BA0400 HORIZONTAL SIDE RAILS (SET)
- ✓ BA0200 200MM CLIP-ON BASE EXTENSION
- ✓ FREE STANDING IV POLE
- ✓ BA0580 CLAMP ON RAIL



Self help pole included  
Full function handset with lockout function

