



www.jbmedical.com.au

JB MEDICAL EQUIPMENT

**DECEMBER 2022**

**BROCHURE**

## ESSEX BEDSIDE COMMODE

**SKU: VR161G**

The Essex Commode Chair is height adjustable, practical commode chair suitable for use as a bedroom commode or bathroom commode. The Essex Commode Chair has a 5 litre bucket with lid and a padded over seat is also supplied for discretion not in use. The upholstery of the commode is well padded and is vinyl welded to reduce the possibility of bacteria build up. Available fully assembled or flat packed.



**\$120**  
ex gst

### Specifications Dimensions

(L x W x H mm):

460 x 540 x 820-1030mm

Seat Length:

420mm

Seat Width:

420mm

Seat Height:

490-700mm

Armrest Height:

170mm

Unit Weight (kg):

6.3kg

Weight Capacity (kg):

158kg

### KEY FEATURES

- PADDED BACK
- DETACHABLE PADDING OVER SEAT
- HEIGHT ADJUSTABLE
- 5 LITRE PAIL WITH LID SUPPLIED
- HINGED THERMOSET PLASTIC TOILET SEAT
- CAN BE DISMANTLED FOR TRANSPORTATION OR STORAGE
- WEIGHT CAPACITY 158KG

**aidapt**

**ON SALE**  
ORDER NOW



I ♥  
ndis

NDIS REGISTERED  
SERVICE PROVIDER



**1300 788 182**