



www.jbmedical.com.au

JANUARY 2023

BROCHURE

WIRELESS BED SENSOR MAT

BED MAT & TRANSMITTER: MED-B365X-WT

FULL KIT: RECEIVER UNIT, BED MAT & TRANSMITTER: MED-B365K-W6T(X)ORW3T(X)

Indigo Care provides a New & Improved pressure sensitive Premium Bed Sensor Mat that is more comfortable & less noisy for the user. The mat's thickness has increased to 12mm providing a softer touch and feel for the user at the same time ensuring the noise created due to contact is quieter. This improved Bed Mat will be a huge hit with your residents and staff! MORE COMFORT, LESS NOISE Our Bed Mats are designed to immediately alert staff should a resident leave their bed which is most useful to residents who are at risk of falling out of their bed, or those who require assistance in arising from the bed. The transmitter is an exclusive ID encoded RF Transmitter which communicates with the Wireless Receiver Box (RRU). It allows the bed sensor mat to wirelessly transmit an alarm to the call point. Simply connect the Transmitter to the bed sensor mat, then teach in the Transmitter to RRU by following the instructions printed. The sensor mat will now trigger the call point when activated within the transmitting area. It can be integrated to most nurse call systems. (Note: Transmitter and Bed mat can be ordered separately)

**MORE COMFORT,
LESS NOISE**



KEY SPECIFICATIONS:

Bed mat

- Dimensions: 760mm X 510mm
- Connection: RJ11 on 500mm curly cord
- Electrical resistance < 500 Ohm
- White PVC, Thickness 12mm
- Warranty: 12 months

Transmitter

- Connection: RJ12 socket
- No External Antenna
- LED battery indicator
- Transmits up to 30m (line of sight)
- Frequency: 433.92MHz
- Uses a 3V battery-CR2032
- Warranty: 12 months

KEY FEATURES

- ✓ MORE COMFORT LESS NOISE
- ✓ INTEGRATES TO MOST NURSE CALL SYSTEMS VIA OUR RECEIVERS
- ✓ NO WIRES AROUND RESIDENT/PATIENT NO TRIP HAZARDS
- ✓ REDUCE THE RISK OF FALLS & INJURIES
- ✓ SOLVENT WIPE-DISINFECTANT
- ✓ LOW PROFILE, NONSLIP MATERIAL
- ✓ PORTABLE
- ✓ EASY TO MOVE FROM ROOM TO ROOM
- ✓ EASY TO MOVE FROM ROOM TO ROOM
- ✓ PLUG & PLAY OPERATION
- ✓ DURABLE & RELIABLE
- ✓ COST EFFECTIVE



**ON SALE
ORDER NOW**

WIRELESS BED SENSOR PAD (UNDER MATTRESS)

BED MAT & TRANSMITTER: MED-B365X-UM-WT

FULL KIT: RECEIVER UNIT, BED MAT & TRANSMITTER: MED-B365K-UM-WT(X)

Indigo Care provides a pressure sensitive Wireless Bed Sensor Pad that is designed for UNDER MATTRESS to immediately alert staff when the resident leaves their bed. Ideally suited for extreme falls risk residents when placed under the shoulder area, it will notify through your nurse call system when the resident is awake and sits up in bed. The transmitter is an exclusive ID encoded RF Transmitter which communicates with the Wireless Receiver Box (RRU) or the Universal Call Point (UCPV5.1). It allows the bed sensor mat to wirelessly transmit an alarm to the call point. Simply connect the Transmitter to the bed sensor mat. Then teach in the Transmitter to RRU or Dual Function Call Point by following the instructions printed. The sensor mat will now trigger the call point when activated within the transmitting area. (Note: Transmitter and Bed mat can be ordered separately) Manufacturer recommendation: Check the sensor mat to see if compatible with the Under Mattress Sensor Pad, before purchasing.



KEY FEATURES

- ✓ SUITABLE FOR UNDER MATTRESS
- ✓ REDUCE THE RISK OF FALLS & INJURIES
- ✓ INTEGRATES TO MOST NURSE CALL SYSTEMS
- ✓ EASY TO MOVE FROM ROOM TO ROOM
- ✓ SOLVENT WIPE-DISINFECTANT
- ✓ DURABLE & RELIABLE
- ✓ LOW PROFILE, NON-SLIP MATERIAL
- ✓ PORTABLE
- ✓ COST EFFECTIVE
- ✓ EASY TO USE
- ✓ PLUG & PLAY OPERATION
- ✓ LONG PRESSURE MAT (760MM X 110MM)

KEY SPECIFICATIONS

Bed Mat

- Dimensions: 760mm X 110mm
- Electrical resistance < 100 Ohm
- White PVC, Thickness 0.20mm
- Warranty: 12 months

Transmitter

- Connection: RJ12 socket
- No External Antenna
- LED low battery indicator
- Transmits up to 30m (line of sight) Frequency: 433.92MHz
- Uses a 3V battery-CR2032 which will last up to 12 months
- Warranty: 12 months

PRODUCT ACCESSORIES

- Wireless Receiver (RRU)
- Wireless Floor Sensor Mat
- Wireless Bed Sensor Mat
- Wireless Under Mattress Mat
- Wireless Chair Mat
- Wireless Waterproof Neck or Wrist Pendant
- Wireless PIR Motion Sensor



NDIS REGISTERED SERVICE PROVIDER



1300 788 182



www.jbmedical.com.au

JANUARY 2023**BROCHURE****WIRELESS CHAIR MAT****CHAIR MAT & TRANSMITTER: MED-C365X-WT****FULL KIT: RECEIVER UNIT, CHAIR MAT & TRANSMITTER: MED-C365K-WT(X)**

Indigo Care provides a pressure sensitive Wireless Chair Sensor Pad that is designed to immediately alert staff when a resident leaves their chair. The transmitter is an exclusive ID encoded RF Transmitter which communicates with the Wireless Receiver Box (RRU). It allows the Chair sensor mat to wirelessly transmit an alarm to the call point. Simply connect the Transmitter to the chair sensor mat. Then program the Transmitter to the RRU by following the instructions printed. The sensor mat will now trigger the call point when activated within the transmitting area. It is easily integrated to most nurse call systems. (Note: Transmitter and Bed mat can be ordered separately)

EASY TO USE**KEY SPECIFICATIONS:****Chair mat**

- Dimensions: 380mm X 255mm
- Connection: RJ12 on 500mm curly cord
- Silver Printing, Resistance <1k Ohm • White PVC, Thickness 0.20mm
- Warranty 12 months

Transmitter

- Connection: RJ 4P4C socket
- No External Antenna
- LED low battery indicator
- Transmits up to 30m (line of sight)
- Frequency: 433.92MHz
- Uses a 3V battery-CR2032
- Warranty: 12 months

PRODUCT ACCESSORIES:

- Wireless Receiver (RRU)
- Wireless Floor Sensor Mat
- Wireless Bed Mat
- Wireless Under Mattress Mat
- Wireless Waterproof Neck Pendant
- Wireless Waterproof Wrist Pendant Wireless PIR Motion Sensor

KEY FEATURES

- ✓ INTEGRATES TO MOST NURSE CALL SYSTEMS VIA OUR RECEIVERS
- ✓ REDUCE THE RISK OF FALLS & INJURIES
- ✓ SOLVENT WIPE-DISINFECTANT
- ✓ LOW PROFILE, NONSLIP MATERIAL
- ✓ NO WIRES AROUND RESIDENT/PATIENT NO TRIP HAZARDS
- ✓ PORTABLE

- ✓ EASY TO USE
- ✓ EASY TO MOVE FROM ROOM TO ROOM
- ✓ PLUG & PLAY OPERATION
- ✓ DURABLE & RELIABLE
- ✓ COST EFFECTIVE

**ON SALE**
ORDER NOWNDIS REGISTERED
SERVICE PROVIDER**1300 788 182**